



## UKST130SS - UKh2o 13 Litre Quick Release Stainless Steel Sediment Trap Kit



The UKh2o UKST130SS all stainless steel (EN 1.4301 grade) sediment traps are constructed for durability, hygiene and professional appearance. Suitable for laboratories, hospitals, dental labs etc. These traps simply screw straight on to the Ø38mm sink drain spigot - what could be simpler?

- A quick-release container allows quick and easy disposal of captured materials
  - Ø38mm threaded top inlet and reversible left or right outlet fittings
  - 13 litre sediment capacity
  - A unique combined release-bar and handle ensures rapid, safe waste disposal
- UKST130SS (13 litre) stainless steel sediment trap provides highly effective separation of finer residues but also heavy metals. Additional wall support is provided via adjustable cantilever brackets.

We recommend emptying the bucket when it is filled to  $\frac{1}{3}$  of its capacity. When emptying, open the tap at the front and drain the water, then close the tap when the water has stopped running. The bucket is released by lifting the handles upwards. The lid will remain attached to the wall and the bucket will drop about 20mm. Slide the bucket towards you and use the handle to carry it. Empty the material from the bucket and reassemble the trap ensuring that the runners touch the stops at the end of the rails.

Installation Kit includes:

- 1 ½ " Universal Outlet Adaptor
- Full Fitting Instructions

• Tel: 01903 500551

• E-Mail: [enquiries@ukh2o.co.uk](mailto:enquiries@ukh2o.co.uk) or visit our web site: [www.ukh2o.co.uk](http://www.ukh2o.co.uk)